


Case Overview

The case details below are courtesy of Dr. Markus Möhlenbruch and his team at University Hospital Heidelberg in Germany*

OPHTHALMIC ICA ANEURYSM

A 63-year-old female presented with a right ophthalmic ICA aneurysm.




Pre-procedural images and measurements are shown above. A FRED X 4 x 18/12mm device was chosen to ensure ample wall apposition (~0.5mm greater than the largest vessel segment).


INTRA-PROCEDURAL DETAILS

Dr. Möhlenbruch used an access setup of Sofia 6F and Headway 27 catheters.

The images below show the FRED X device immediately post deployment. Dr. Möhlenbruch used a two-handed deployment technique where he predominantly unsheathed the microcatheter while holding steady forward pressure on the device delivery wire. Dr. Möhlenbruch was careful to keep the system centered in the vessel throughout the deployment.



Immediately post deployment, the device looked completely open and well apposed to the vessel wall. Dr. Möhlenbruch removed the delivery wire slowly and carefully to ensure that the proximal finished ends did not get caught on the distal tip of the delivery wire.




There were no DWI lesions present immediately post deployment.

*THIS IS A CASE STUDY EXAMPLE. INDIVIDUAL PATIENT RESULTS MAY VARY.

DUAL ANTIPLATELET (DAPT) REGIMEN


Dr. Möhlenbruch prescribed a DAPT regimen of Aspirin (on-going) and Clopidogrel (for 3 months).

3 MONTH FOLLOW-UP




At three-month follow-up, the aneurysm is partially occluded and the device remains fully open, with no complications.

6 MONTH FOLLOW-UP




At six-month follow-up, there is complete aneurysm occlusion.

12 MONTH FOLLOW-UP



Twelve-month follow-up continues to show complete aneurysm occlusion.

24 MONTH FOLLOW-UP



At twenty-four-month follow-up, results are consistent and the aneurysm continues to be completely occluded.

***THIS IS A CASE STUDY EXAMPLE. INDIVIDUAL PATIENT RESULTS MAY VARY.**



For more information, contact your local MicroVention sales representative or visit our website. www.fred-x.com

FRED™ X System Product Indications, U.S.: The Flow Re-Direction Endoluminal Device (FRED) System is indicated for use in the internal carotid artery from the petrous segment to the terminus for the endovascular treatment of adult patients (22 years of age or older) with wide-necked (neck width \geq 4 mm or dome-to-neck ratio $<$ 2) saccular or fusiform intracranial aneurysms arising from a parent vessel with a diameter \geq 2.0 mm and \leq 5.0 mm.

FRED™ X System Product Indications, EMEA: The FRED system is intended for endovascular embolization of intracranial neurovascular aneurysms. The FRED system may also be used with embolic coils for the treatment of intracranial neurovascular lesions.

HEADWAY™ Microcatheter Indications for Use: The HEADWAY Microcatheter is intended for general intravascular use, including the peripheral, coronary and neuro vasculature for the infusion of diagnostic agents, such as contrast media, and therapeutic agents, such as occlusion coils.

SOFIA™ EX Catheter Indications for Use: The SOFIA EX Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA EX Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA EX Catheter is not intended for use in coronary arteries.

SOFIA™ Catheter Indications for Use: The SOFIA Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA Catheter can be used to facilitate introduction of diagnostic or therapeutic agents. The SOFIA Catheter is not intended for use in coronary arteries. Moreover, the SOFIA Catheter is intended for use in removal/aspiration of emboli and thrombi from selected blood vessels in the arterial system, including the peripheral and neuro vasculatures.

For Healthcare Professional Use Only

RX Only: Federal (FDA) law restricts this device to sale by or on the order of a physician.

MicroVention, FRED, stylized X, SOFIA, and HEADWAY are registered trademarks of MicroVention, Inc. in the United States and other jurisdictions. Third party brands are trademarks of their respective owners. For complete indications, contraindications, potential complications, warnings, precautions, and instructions, see instructions for use (IFU provided in the device).

©2023 MicroVention, Inc. MM1495(j) EMEA-US 05/23

CE0297 class III

**MicroVention Worldwide
Innovation Center** 1.714.247.8000
35 Enterprise
Aliso Viejo, CA 92656 USA
Customer Service 1.800.990.8368
Website microvention.com

MicroVention UK Limited
Suite 3, The Barrack Building
10 Cliffords Fort, North Shields
Tyne and Wear, NE30 1JE
United Kingdom
PH +44(0) 191 258 6777
F +44(0) 191 258 5999

MicroVention Europe, S.A.R.L.
(Legal Manufacturer EMEA)
30 bis, rue du Vieil Abrevoir
78100 Saint-Germain-en-Laye
France
RCS 440 775 674
with capital of 40 000€
PH +33(1) 39 21 77 46
F +33(1) 39 21 16 01

MicroVention Deutschland GmbH
Hildebrandtstr. 4 F
D-40215 Dusseldorf
Germany
PH +49 211 210 798-0
F +49 211 210 798-29

MicroVention Italia S.r.l.
Via Tommaso Gulli n. 39
20147 Milano
Italy
PH +39 071 7106156
F +39 071 2865400

